



# Saving Lives Across Australia

with

## Heart Safe Hotel

Are You

**Ready, Willing & Able to HELP?**



## Sudden Cardiac Arrest (SCA) knows NO Boundaries

In Australia, approximately 30,000 people each year sustain SCA out-of-hospital and are treated by emergency medical services (EMS).

**That's 575 Australians suffering SCA every week... and ONLY 1 out of 10 WILL SURVIVE!**



### The Facts

- ♥ On average, only 10% of SCA victims in Australia survive
- ♥ It can happen **anytime, anywhere and to anyone...** even fit, young athletes
- ♥ **Early defibrillation is the single most effective treatment for SCA**
- ♥ **Defibrillation within three minutes of collapse can increase the chance of survival to over 70%**



Raising Awareness & Education  
About Saving Lives  
in Your Hotel

## Three Steps to Save a Life



**1. CALL 000**



**2. PUSH**



**3. SHOCK**



**\$1,795**

(tax deductible item)

[\*\*ORDER NOW\*\*](#)



**\$1,995 (Bundle)**

(tax deductible item)

[\*\*ORDER NOW\*\*](#)

**Mindray C1 Public Defibrillator designed to address challenges faced by first-time rescuers by properly guiding users through the rescue process**

- ♥ Only 1 set of pads with a switch for **Paediatric Mode**
- ♥ Provides escalating energy up **from 200J up to 360J**
- ♥ **4+-Year** life cycle for pads & battery
- ♥ **Fastest time to first shock** (< 8 seconds)
- ♥ 1.5-meter drop endurance
- ♥ Highest rated for Dust and Water Protection (IP55)

**Adheres to Australian Resuscitation  
Council Guidelines (ANZCOR)**

Contact [Heart Safe Australia](#) to discuss any questions about a Heart Safe Hotel: